# Rosa de Bengala

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1. Young E J. Clinical Infectious Diseases 1995; 21: 283-290.
2. Alton GC. Techniques for Brucellosis Laboratory INRA Paris, 1988.
3. Ariza J. Current Opinion in Infectious Diseases 1996; 9: 126-131.
4. Comité mixto FAO/OMS de expertos en Brucelosis. WLD Health Org Tech Rep Ser 1958; 148: 1-60.
5. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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|              |          |                            |
|--------------|----------|----------------------------|
| Ref. 1200901 | 50 tests | : 2,5 mL Rosa de Bengala   |
|              |          | : 1 mL Control +           |
|              |          | : 1 mL Control -           |
|              |          | : 8 x 6 portas desechables |

**Qualitative determination of antibodies anti-Brucella IVD**

Store at 2 - 8°C.

**PRINCIPLE OF THE METHOD**

The Rose Bengal is a slide agglutination test for the qualitative and semi-quantitative detection of antibodies anti-Brucella in human serum. The stained bacterial suspension agglutinates when mixed with samples containing specific IgG or IgM antibodies present in the patient sample.

**CLINICAL SIGNIFICANCE**

Brucella diagnostic may be assessed either by microorganism isolation in blood or stools, or by titration of specific antibodies in the patient serum. The reagent, because of its formulation in an acid buffer, is reactive with both IgG and IgM antibodies and very useful for the diagnosis of chronic individuals, which present a high level of IgG antibody, difficult to be detected by the reference tube method (Wright).

**REAGENTS**

|                                 |  |
|---------------------------------|--|
| <b>Rose Bengal</b><br>White cap | <i>Brucella abortus</i> suspension, strain S99, in lactate buffer 1 mol/L, phenol 5 g/L, Rose Bengal, pH 3,6 |
| <b>Control +</b><br>Red cap     | Animal serum, with an antibody anti- <i>Br.abortus</i> concentration $\geq$ 50 IU/mL. Preservative           |
| <b>Control -</b><br>Blue cap    | Animal serum. Preservative   |

**CALIBRATION**

The Rose Bengal sensitivity is calibrated against the 2<sup>o</sup> International Preparation of anti-*Brucella abortus* from NIBS (UK) (WHO).

**STORAGE AND STABILITY**

All reagents are ready to use, and will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze: frozen reagents could change the functionality of the test. Always keep vials in vertical position. If the position is changed, gently mix to dissolve aggregates that may be present.

**Reagents deterioration:** Presence of particles.

**ADDITIONAL EQUIPMENT**

- Mechanical rotator with adjustable speed at 80-100 r.p.m.
- Vortex mixer.
- Pippetes 50  $\mu$ L.

**SAMPLES**

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before use. Do not use highly hemolized or lipemic samples.

**PROCEDURE**
**Qualitative method**

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50  $\mu$ L of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
3. Mix the R. Bengal reagent vigorously or on a vortex mixer before using and add one drop next to the sample to be tested.
4. Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
5. Place the slide on a mechanical rotator at 80-100 r.p.m. for 4 minutes. False positive results could appear if the test is read later than two minutes.

**Semi-quantitative method**

1. Make serial two fold dilutions of the sample in 9 g/L saline solution.
2. Proceed for each dilution as in the qualitative method.

**READING AND INTERPRETATION**

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator. The presence of agglutination indicates an antibody anti-Brucella concentration equal or greater than 25 IU/mL. The titer, in the semi-quantitative method, is defined as the highest dilution showing a positive result.

**CALCULATIONS**

The approximate antibody concentration in the patient sample is calculated as follows:

$$25 \times \text{anti-Brucella Titer} = \text{IU/mL}$$

**QUALITY CONTROL**

Positive and Negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation. All result different from the negative control result, will be considered as a positive.

**REFERENCE VALUES**

Up to 25 IU/mL.

Each laboratory should establish its own reference range.

**PERFORMANCE CHARACTERISTICS**

1. **Analytical sensitivity:** 25 ( $\pm$ 5) IU/mL, under the described assay conditions
2. **Prozone effect:** No prozone effect was detected up to 1000 IU/mL.
3. **Diagnostic sensitivity:** 100 %.
4. **Diagnostic specificity:** 98 %.

**INTERFERENCES**

Hemoglobin (10 g/L), lipemia (10 g/L) and rheumatoid factors (300 IU/mL), do not interfere. Bilirubin interferes at 2.5 mg/dL. Other substances may interfere<sup>5</sup>.

**NOTES**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**BIBLIOGRAPHY**

1. Young E J. Clinical Infectious Diseases 1995; 21: 283-290.
2. Alton GC. Techniques for Brucellosis Laboratory INRA Paris, 1988.
3. Ariza J. Current Opinion in Infectious Diseases 1996; 9: 126-131.
4. Comité mixto FAO/OMS de expertos en Brucelosis. WLD Health Org Tech Rep Ser 1958; 148: 1-60.
5. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995

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