BILIRUBIN D- DMSO

Bilirubin Direct

DMSO. Colorimetric

Quantitative determination of bilirubin IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Bilirubin is converted to colored azobilirubin by diazotized sulfanilic acid and measured photometrically. Of the two fractions presents in serum, bilirubin-gluconamide and free bilirubin loosely bound to albumin, only the former reacts directly in aqueous solution (bilirubin direct), while free bilirubin requires solubilization with dimethylsulphoxide (DMSO) to react (bilirubin indirect). In the determination of indirect bilirubin the direct is also determined; the results correspond to total bilirubin. The intensity of the color formed is proportional to the bilirubin concentration in the sample^{1,2,3}.

CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of hemoglobin. It is transported from the spleen to the liver and excreted into bile. Hyperbilirubinemia results from the increase of bilirubin concentrations in plasma.

Causes of hyperbilirubinemia:

Total bilirubin: Increase hemolysis, genetic errors, neonatal jaundice, ineffective erythropoiesis, and drugs.

Direct bilirubin: Hepatic cholestasis, genetic errors, hepatocellular damage^{1,6,7}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Sulfanilic acid	30 mmol/L
	Hydrochloric acid (HCl)	150 mmol/L
R 2	Sodium nitrite	29 mmol/L

PRECAUTIONS

R1: H290-May be corrosive to metals. H314-Causes severe burns and eye damage. EUH208-Contains sulphanic acid. May produce an allergic reaction.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Pipette 1,5 mL of R2 into R1 content. Mix avoiding foam forming and it will be ready to use (WR). Do not use the reagent before 30 min. after the reagent preparation.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Color development in R 2.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyser.
- General laboratory equipment

SAMPLES

Serum or plasma, free of hemolysis¹. Protect samples from direct light. Stability: Bilirubin is stable at 2-8°C for 4 days and 2 months at -20°C.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	250
Test	BILD	Vol. R2	
Full Name	Direct Bilirub.	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	20	Decimal	0.01 Slope 1
Vol. Sample Increas.		Unit	mg/dL Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	546	Reagent Blank	10-11
Sec. Wave.		React. Time	46-47

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	BILD	No.	**
Full name	BILD	Print name	BILD
Reac. Type	End Point	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	mg/dL	Decimal	0.01
Reagent Blank	10 - 11	React. Time	46 - 47
Vol. Sample	20 ul	R1	250 ul
Increased		R2	
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **28 days**. After this period or in case of bad results, the working reagent must be prepared again.

REFERENCE VALUES¹

Bilirubin Direct: Up to 0,25 mg/dL \approx 4,27 μ mol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,07 mg/dL to linearity limit of 20 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	0,96	2,48	0,96
SD	0,024	0,051	0,043	0,035
CV (%)	2,52	2,06	4,49	1,41

Sensitivity: 1 mg/dL = 0,006856 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,96.

Regression equation: y=0,71177x - 0,05267

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
2. Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
3. Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAPHIE

1. Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
2. Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
3. Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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