



CE BILIRUBIN T-DMSO

Bilirubin Total

DMSO. Colorimetric

Quantitative determination of bilirubin

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Bilirubin is converted to colored azobilirubin by diazotized sulfanilic acid and measured photometrically. Of the two fractions present in serum, bilirubin-glucuronide and free bilirubin loosely bound to albumin, only the former reacts directly in aqueous solution (pH=6.8), while the latter requires solubilization with dimethylsulfoxide (DMSO) to react (bilirubin indirect). In the determination of indirect bilirubin the direct is also determined, the results correspond to total bilirubin. The intensity of the color formed is proportional to the bilirubin concentration in the sample.^{1,2,3}

CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of hemoglobin. It is transported from the spleen to the liver and excreted into bile. Hyperbilirubinemia results from the increase of bilirubin concentrations in plasma. Causes of hyperbilirubinemia: Total bilirubin: increase hemolysis, genetic errors, neonatal jaundice, ineffective erythropoiesis, and drugs. Direct bilirubin: Hepatic (cholestasis, genetic errors, hepatocellular damage)^{1,2}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

	SUBSTRATE	30 mmol/L
R 1	diazotized sulfanilic acid (HCl)	50 mmol/L
	dimethylsulfoxide (DMSO)	7 mmol/L
R 2	Sodium nitrite	30 mmol/L
Optional	BILIRUBIN CAL	Ref. 100250

PRECAUTIONS

R1/R2: Corrosive (C) R55: Causes severe burns. S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contamination prevention during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Color development in R2.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 555 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis¹. Protect samples from direct light. Stability: Bilirubin is stable at 2-8°C for 4 days and 2 months at -20°C.

PROCEDURE

1. Assay conditions:

Wavelength: 555 nm (535-580)
Cuvette: 1 cm light path
Temperature: 15-25°C

2. Adjust the instrument to zero with distilled water.

3. Pipette into a cuvette:

	Blank	B. Total
R 1 (mL)	1.5	1.5
R 2 (μL)	40	40
Sample ^{1,2,3} / Calibrator (μL)	100	100

4. Mix and incubate for exactly 5 minutes at room temperature.

5. Read the absorbance (A).

CALCULATIONS

- With Calibrator:

$$\frac{(A) \text{ Sample}}{(A) \text{ Calibrator}} \times \frac{(A) \text{ Sample Blank}}{(A) \text{ Calibrator Blank}} \times \text{Conc. Calibrator} = \text{mg/dL bilirubin}$$

- With Factor:

$$(A) \text{ Sample} - (A) \text{ Sample Blank} \times \text{Factor}^* = \text{mg/dL bilirubin in the sample}$$

$$\text{Factor} = \frac{\text{Concentration of Calibrator}}{(A) \text{ Calibrator} - (A) \text{ Calibrator Blank}}; \text{Theoretical factor} = 10.1$$

Conversion factor: mg/dL $\times 17.1 = \mu\text{mol/L}$

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINREACT H Normal and Pathologic (Ref. 1002120 and 1002130).

Control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Bilirubin Total Up to 1.10 mg/dL $\approx$ Up to 18.81 $\mu\text{mol/L}$

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0.099 mg/dL to theory limit of 16 mg/dL.

1. The results obtained were greater than linearity limit, since the sample 1:2 with NaOH pH, and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	1.12	0.36	1.11	0.35
SD	0.05	0.12	0.03	0.12
CV (%)	4.35	32.7	2.70	33.2

Sensitivity: 1 mg/dL $\approx$ 0.01340 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

Correlation coefficient (r): 0.998.

Regression equation: $y = 0.0206 + 0.9884 x$.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Hemolysis causes decreased bilirubin values.^{2,3}

1. All of drugs and other interfering substances with bilirubin has been reported.¹

NOTES

- For bilirubin determination is necessary, pipette 50 μL of sample. Multiply the result by 2.

1. SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Kasari A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St. Louis, Toronto, Princeton 1984; 32:9-124; 418 and 630.
- Melchior H T et al. The determination of bilirubin with the photometric colorimeter. J. Biol Chem 1937; 112: 248-491.
- Machuga, R. Improved micro-method for determination of serum bilirubin. Clin Chem 1966; 12: 15-17.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed M&C Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed M&C Press, 1995.
- Burke A et al. Textbook of Clinical Chemistry, 3rd ed AACO 1988.
- Tietz NW et al. Clinical Guidelines Laboratory Tests, 3rd ed AACO 1995.

PACKAGING

Ref: 1001042
R 1: 2 x 150 mL
R 2: 1 x 10 mL

BSS04-I 02/03/2

SPINREACT S.A.S. - Chem. Santa Colomba, 7 E-17120 SANTA ESTEVE DE BAS (GI) SPAIN
Tel: +34 932 93 06 80 Fax: +34 932 93 00 19 e-mail: spinfo@spinfo.com